



PRODUCT SPECIFICATION



Product: CEDERROTH EYE&WOUND CLEANSING
SPRAY 150 ML

Reference code: 726000



1. General description

A portable spray bottle filled with 150 mL physiological (0.9%) unbuffered sodium chloride solution, preservative free. The nozzle is suitable for use on both eyes and wounds. The actuator is protected by a plastic cap and a tamper seal. The nitrogen propellant is not in contact with the solution. The solution remains sterile until last drop.

2. Intended use

The Cederroth Eye&Wound Cleansing spray is intended to be used for rinsing and cleansing out dust and dirt from eyes and wounds.

3. Technical description

Product dimensions

Spray bottle (can)

Width (mm): 53
Height (mm) 162
Fill volume (ml): 150

Material

Liquid Sterile 0.9% Sodium Chloride (NaCl). Preservative free.
Can Aluminum.
Bag-on-valve 4 layer PP-OPA-ALU-PET
Actuator White POM
Cover CAP Transparent PP

4. Sterilization

Sterilization method: Irradiation

5. Composition

Consumer package

<u>Ref#</u>	<u>Product/Material</u>	<u>Pcs</u>
726000	Cederroth Eye&Wound Cleansing Spray	1

Transport package

Number of consumer packages in transport package:	12
Dimensions (l x w x h) mm:	221 x 168 x 191
Weight (kg):	2,5



PRODUCT SPECIFICATION



6. Labeling

On product = consumer package

Cederroth contact information
Product name, Volume, Reference code,
CE+NB#, SterileR, Expiry date, LOT, EAN-code,
Temperature limitations, Instructions for use, 3

On transport package

Cederroth contact information
Product name, Reference code
Lot-number, Expiry date, EAN-
code, Limited quantity symbol
Quantity, CE+NB#, SterileR

7. Storage conditions

Store in temperatures between 0°C and +40°C.
Shelf life: 5 years from production date

Authorization: